

Ingenuity in  
infection control  
since 1830



# Instrument Management Systems for Eschmann Autoclaves & Washer Disinfectors.

**Instrument  
Management Systems**

**Exclusive  
to Eschmann**

# Instrument Management Systems



**Eschmann has a range of instrument management systems that are simple, effective and most importantly compliant, helping to suit the needs of your clinical environment.**

The range is available in different colours to help practitioners easily identify their instruments and adopt instrument management traceability in accordance with guidance.

## Instrument Trays

**Perforated aluminium and stainless steel trays, suitable for use in vacuum and non-vacuum autoclaves are available in three sizes - Standard (28cm x 18cm), Examination (18cm x 14cm) and a long tray specifically designed for use with the SES 3020B autoclave (38cm x 19cm).**

Examination trays hold enough instruments for regular check-up procedures and can be assigned to individual patients, without compromising the capacity of the autoclave.

Both examination and standard trays are available in a choice of six colours.



### Ordering Information

**Note:** The product code last digit 'x' determines the colour. Replace 'x' with either: 1-Silver, 2- Blue, 3-Red, 4-Green, 5-Yellow, 6-Turquoise

| Part No. | Description  |
|----------|--|
| 872445x  | Standard Tray (28cm x 18cm) Aluminium                  |
| 872446x  | Examination Tray (14cm x 18cm) Aluminium               |
| 8725450  | Standard Tray (28cm x 18cm) Stainless Steel            |
| 8725460  | Examination Tray (14cm x 18cm) Stainless Steel         |
| 8721451  | Standard Storage Lid (28cm x 18cm) Aluminium Silver    |
| 8721461  | Examination Storage Lid (14cm x 18cm) Aluminium Silver |
| 8725470  | Long Tray for SES 3020B (38cm x 19cm) Stainless Steel  |
| 8727121  | Tray Lifter  |

# HFIT Plus Cassettes

**Hands Free Instrument Transfer (HFIT) Plus Cassettes significantly reduce the risk of sharps injury whilst ensuring total cleaning and sterilising efficacy.**

Available in three sizes, these innovative stainless steel, self locking, HFIT Plus Cassettes securely hold instruments enabling them to be safely passed through the complete decontamination and sterilisation process without the need for user handling.

The HFIT Plus Cassette silicones are available in a choice of six colours.



## Ordering Information

**Note:** The product code last digit 'x' determines the colour. Replace 'x' with either: 2- Blue, 3-Red, 4-Green, 5-Yellow, 6-Silver or 8-White

| Part No. | Description           |
|----------|-----------------------|
| 872280x  | HFIT-5 Plus Cassette  |
| 872281x  | HFIT-10 Plus Cassette |
| 872282x  | HFIT-20 Plus Cassette |

# Burr Holders

**Stainless steel burr holders with silicone plugs securely and safely hold any combination of FG or CA burs, reducing the risk of sharps injury during transfer.**

The cover can be used to offer the burs flat or at an angle.

Available in a choice of four colours.



## Ordering Information

**Note:** The product code last digit 'x' determines the colour. Replace 'x' with either: 2-Blue (dental burr holder only), 3-Red, 4-Green, 5-Yellow

| Part No. | Description          |
|----------|----------------------|
| 872021x  | Dental Burr Holder   |
| 872600x  | Surgical Burr Holder |



# Purchase your Eschmann consumables & accessories online at [eschmann.co.uk](https://eschmann.co.uk).

SAVE 10% on all Instrument Management products purchased from the Eschmann website using code:

**ESCHMANN\_IMS10**



**Eschmann Technologies Limited**  
Eschmann House, 15 Peter Road, Lancing,  
West Sussex, BN15 8TJ United Kingdom

**t:** +44 (0)1903 753322  
**e:** [info@eschmann.co.uk](mailto:info@eschmann.co.uk)  
**w:** [eschmann.co.uk](https://eschmann.co.uk)

Ingenuity in  
infection control  
since 1830

